

S series

行号	产品代码	产品名称	产品名称 (英文)	规格型号	单位	图片	S30S	S30X	S60S	S60X	S90S	S90X
1	1.00.120DL002C0	LED手术灯V1	LED Operating light V1 (without /trademark plate)	2P Connector	pcs		✓	✓	✓	✓	✓	✓
2	1.00.120DS002A0	LED手术灯V2(不带透明罩)	LED Operating light V2 (without transparent cover)	2P Connector	pcs		✓	✓	✓	✓	✓	✓
3	1.00.120DB002A2	LED手术灯臂(西格)	LED Operating light arm (Siger)	2P Connector	pcs		✓	✓	✓	✓	✓	✓
4	1.01.111JT001F0	S90踏板脚踏总成	Pedal foot control assembly		pcs						✓	✓
5	1.01.111ZC010F1	90电控脚踏总成(瓷白注塑)	Electric foot control assembly(white,plastic injection)		pcs						✓	✓
6	1.01.112JT001F5	S30踏板脚踏总成	Pedal foot control assembly		pcs		✓	✓	✓	✓		
7	1.01.112JT002F5	S30调速踏板脚踏总成	Speed-regulating pedal foot control assembly		pcs		✓	✓	✓	✓		
8	1.01.112ZC012F5	30/60气控脚踏总成(瓷白注塑)	Pneumatic foot control assembly(white,plastic injection)		pcs		✓	✓	✓	✓		
9	1.01.112ZC013F5	30/60调速脚踏总成(瓷白注塑)	Speed-regulating pedal foot control assembly (white,plastic injection)		pcs		✓	✓	✓	✓		
10	1.01.111ZC001F1	90主操作触摸面板总成(瓷白)	S90 Main control touch panel assembly (white)		pcs						✓	
11	1.01.112ZC001F5	60主操作触摸面板总成(瓷白带显示屏)	Main control touch panel assembly (white, with display screen)		pcs		✓		✓			
12	1.01.112ZC002F5	60助手位触摸上盖总成(瓷白)	Assistant touch top cover assembly (white)		pcs		✓	✓	✓	✓	✓	✓
13	1.01.112ZL003F5	30/60下挂上罩壳总成(瓷白带显示屏)	X-type instrument table top cover assembly (white, with display screen)		pcs			✓		✓		
14	1.01.111ZL010F1	90下挂上罩壳总成(瓷白)	X-type instrument table top cover assembly(white)		pcs							✓
15	3.01.301CX021A2	90助手位触摸面板(A2版)	Assistant touch panel (A2)	11 button	块		✓	✓	✓	✓	✓	✓
16	3.02.302CX001A3	60触摸面板	Touch panel		pcs		✓		✓			
17	3.01.301CX003A2	90主操作触摸面板(A2版)	Main control touch panel (A2)	22 button	pcs						✓	

S series

18	3.02.303CX006A2	30下挂大触摸面板(西格)	X-type big touch panel (Siger)		pcs			✓		✓				
19	3.01.301MB002A0	90下挂大触摸面板	X-type big touch panel	23 button	pcs									✓
20	1.01.111ZC008F3	90平衡臂总成(瓷白)	S90 Balance arm assembly (white)		pcs								✓	✓
21	1.02.117ZC005F2	U200平衡臂总成(瓷白上拉)	S-type Balance arm assembly (white)		pcs			✓	✓	✓	✓			
22	1.01.111ZL004A0	上拉三用枪管总成	S-type syringe tubing assembly		pcs			✓		✓		✓		
23	1.01.111ZL005A0	下挂三用枪管总成	X-type syringe tubing assembly		pcs				✓		✓			✓
24	1.01.111ZL006A0	上拉四孔手机管总成	S-type 4-hole turbine tubing assembly		pcs			✓		✓		✓		
25	1.08.191ZL007A0	上拉两孔手机管总成	S-type 2-hole turbine tubing assembly		pcs			✓		✓		✓		
26	1.01.111ZL007A0	下挂四孔手机管总成	X-type 4-hole turbine tubing assembly		pcs				✓		✓			✓
27	1.08.191ZL008A0	下挂两孔手机管总成	X-type 2-hole turbine tubing assembly		pcs				✓		✓			✓
28	1.01.111ZC022A0	助手位三用枪管总成	Assistant syringe assembly		pcs			✓	✓	✓	✓	✓	✓	✓
29	3.08.301GQ0001	国产光纤线	Domestic F.O. turbine tubing	160*180*190CM	pcs			✓	✓	✓	✓	✓	✓	✓
30	1.07.111ZC022A0	助手位三用枪管总成	Assistant syringe assembly	Without angled syringe	pcs			✓	✓	✓	✓	✓	✓	✓
31	3.08.301AB002A0	光导手机线	Bien Air F.O turbine tubing	BIEN NO:1600097-001	pcs			✓	✓	✓	✓	✓	✓	✓
32	3.08.301AB006A0	TKD四孔光纤手机线(带独立回气管)	TKD 4-hole F.O. turbine tubing (with independent air returning pipe)	322.63	pcs			✓	✓	✓	✓	✓	✓	✓
33	3.08.301AG001A0	啄木鸟光固化机(10秒)	Woodpecker curing light (10 seconds)		pcs			✓	✓	✓	✓	✓	✓	✓
34	2.01.201AZ001D3	90器械盘控制板DSP	Instrument table control board (DSP)	90ZL-01/D3	pcs								✓	✓

S series

35	2. 01. 201AZ002D3-1	90器械盘控制主板	Instrument table main control board	90ZL-02/D3(No battery)	pcs								✓	✓
36	4. 01. 418DC20323	钮扣电池	Button battery	LIR2032 Rechargeable battery	pcs								✓	✓
37	2. 01. 203CX001F3A	30侧箱控制板	S30 Hydro group control board	30CX-01-F1	pcs								✓	✓
38	2. 01. 201AZ003D4	90高速手机控制板	Turbine control board	90ZL-03/D4	pcs								✓	✓
39	2. 01. 201AZ005D2	90马达连接板	Motor connection board	90ZL-05/D2	pcs								✓	✓
40	2. 01. 201AZ011D4	90洁牙机控制板	Scaler control board	90ZL-11/D4	pcs								✓	✓
41	2. 01. 201AZ005F1	90马达连接板	Motor connection board	BLDC-05-F1	pcs								✓	✓
42	2. 01. 201AZ006D6	90主操作触摸控制板(大)	Main control touch control board (Big)	90ZL-06/D6	pcs								✓	
43	2. 01. 201AZ016D7	90主操作触摸控制板(小)	Main control touch control board (Small)	90ZL-16/D7	pcs								✓	
44	2. 01. 201AZ008D3	90气刹控制板	Air brake control board	90ZL-08/D1	pcs								✓	
45	2. 01. 201AZ009A2	90新下挂器械盘触摸控制板(大)	X-type Instrument table touch control board (Big)	90ZL-09/A2	pcs									✓
46	2. 01. 201AZ010A3	90新下挂器械盘触摸控制板(小)	X-type Instrument table touch control board (Small)	90ZL-10/A3	pcs									✓
47	2. 01. 201AZ015A2	90新下挂器械盘气刹板	X-type Instrument table air brake(New)	90ZL-15/A0	pcs									✓
48	4. 02. 422LM30333	液晶块	LCD block	LM3033CDW-OB-1 3.3V(ET-12864FV1-GBSW) white with blue word	pcs								✓	✓
49	2. 01. 220DL100G1	LED手术灯主控电路板	LED operating light main control board	LED-01-G1	pcs		✓	✓	✓	✓	✓	✓	✓	✓
50	2. 01. 220DL400G1	LED手术灯电源板	LED Operating light power panel	LED-05-G1	pcs		✓	✓	✓	✓	✓	✓	✓	✓
51	2. 01. 220DL200F5	LED手术灯光源电路板	LED operating light illuminant control board	LED-02-F5(4300K)	pcs		✓	✓	✓	✓	✓	✓	✓	✓
52	2. 01. 220DS003F0	反射式LED光源板	Reflective LED light source board	LED-03-F0	pcs		✓	✓	✓	✓	✓	✓	✓	✓

S series

53	2.01.220DS004F2	反射式LED控制板	Reflective LED control board	LED-04-F2	pcs		✓	✓	✓	✓	✓	✓
54	2.01.201JT002F7	90脚踏方向控制板 (lever type)	Foot pedal control board for chair movement (Lever)	90JT-02-F6	pcs		✓	✓	✓	✓	✓	✓
55	2.01.251JT001F5	V1000脚踏方向控制板 (petal type)	Foot pedal control board for chair movement (Lever)	V1000JT-01-F5	pcs		✓	✓	✓	✓	✓	✓
56	2.01.201AZ004D3	90马达控制板	Motor control board	90ZL-04/D3	pcs		✓	✓	✓	✓	✓	✓
57	2.01.201AZ004F0	90马达控制板 (TKD)	Motor control board (TKD)	BLDC-F0	pcs		✓	✓	✓	✓	✓	✓
58	2.01.201AZ004F0A	MCX马达控制板 (Bien Air)	MCX Motor control board (Bien Air)	BLDC-F0A	pcs		✓	✓	✓	✓	✓	✓
59	2.01.201ZL013F5	90下挂器械盘光控板	X-type instrument table optical control board	90ZL-13-F5	pcs		✓	✓	✓	✓	✓	✓
60	2.01.201ZS001F5	90助手位控制板	Assistant control board	90ZS-01/F5	pcs		✓	✓	✓	✓	✓	✓
61	2.01.201ZS002F5	90助手位触摸控制板	Assistant control board (touch version)	90ZS-02/F5	pcs		✓	✓	✓	✓	✓	✓
62	2.01.203ZL002F5	30新下挂器械盘触摸控制板	New X-type instrument table control board (touch version)	30ZL-02/F5	pcs		✓	✓	✓	✓	✓	✓
63	2.01.203ZL003F5	30主操作触摸控制板	Main control touch control board	30ZL-03/F5	pcs		✓	✓	✓	✓	✓	✓
64	2.01.250GK004F5	U500助手位强吸光控板	Assistant HV suction optical control board	U500ZS-GK-F5	pcs		✓	✓	✓	✓	✓	✓
65	2.01.250GK005F5	U500助手位弱吸光控板	Assistant saliva ejector optical control board	U500ZS-GK-F5	pcs		✓	✓	✓	✓	✓	✓
66	2.01.203CX001F6-A5	30侧箱控制板 (红色J12)	Hydro group control board (Red J12)	30CX-01-F6/A4	pcs		✓	✓	✓	✓	✓	✓
67	2.01.203ZL001F6-A4	30器械盘控制板	Instrument table control board	30ZL-01-F6	pcs		✓	✓	✓	✓	✓	✓
68	2.01.201ZY001F6	座椅控制板 (通用)	Chair control board (universal)	ZY-01-F6	pcs		✓	✓	✓	✓	✓	✓
69	2.01.203JT001F7	30脚踏磁控板	Magnetic control board for foot control	30JT-01-F7	pcs		✓	✓	✓	✓	✓	✓

S series

70	2.01.201ZY021F5	座椅传感器总成	Backrest sensor assembly		pcs			✓	✓	✓	✓	✓	✓
71	2.01.250ZY022F5	U500椅背传感器总成	Backrest sensor assembly		pcs			✓	✓	✓	✓	✓	✓
72	4.02.423AA12060	90升降传动器(不带记忆)	Chair motor without feedbacks	LA31-C215-01 /6000N,120/325,vertical,0.5 meter	pcs			✓	✓	✓	✓	✓	✓
73	4.02.423AA18040	90背部传动器(不带记忆/快速马达)	Backrest motor(without feedbacks/quick motor)	LA31-C123-00 /4000N,180/373,horizontal,1.5 meter	pcs			✓	✓	✓	✓	✓	✓
74	1.01.112ZC003A1	60气减压总成	Air pressure reducing assembly		pcs			✓	✓	✓	✓	✓	✓
75	1.01.112ZC004A1	60气过滤总成	Air filter assembly		pcs			✓	✓	✓	✓		
76	1.01.112ZC005A0	60水过滤总成	Water filter assembly		pcs			✓	✓	✓	✓		
77	1.01.111ZC018A1	90水过滤总成	Water filter assembly		pcs							✓	✓
78	1.01.111ZC019A1	90气过滤总成	Air filter assembly		pcs							✓	✓
79	1.01.111CX020A0	手机回气瓶(加工)	Turbine air returning bottle		pcs			✓	✓	✓	✓	✓	✓
80	1.01.111ZL009F1	拉杆总成(瓷白不带开关)	Pull rod assembly (white, without switch)		pcs			✓		✓		✓	
81	1.01.111ZL008F1	拉杆总成(瓷白带开关)	Pull rod assembly (white, with switch)		pcs			✓		✓		✓	
82	2.03.201ZL214A1-1	90治疗台下盖	Instrument table bottom cover		pcs			✓		✓		✓	
83	2.03.201ZL410A3	90按键面板前后盖	Button panel cover (front & back)		pcs			✓		✓		✓	
84	2.03.201ZL412A0-1	治疗台大面板	Instrument table panel (big)		pcs			✓		✓		✓	
85	2.03.201ZL413A0-1	治疗台小面板	Instrument table panel (small)		pcs			✓		✓		✓	
86	2.06.201ZL426A0	搁架硅胶垫	S-type silicone pad		pcs			✓		✓		✓	












S series

87	2. 03. 201ZL414A0-1	治疗台拉手	Instrument table handle		pcs			✓		✓		✓			
88	2. 03. 201ZL415A0-1	治疗台拉手侧盖	Instrument table handle side cover		pcs			✓		✓		✓			
89	2. 03. 202XZ002A1	大器械盘上罩壳	Big instrument table top cover		pcs				✓		✓		✓		
90	2. 06. 201XZ009A0	大器械盘硅胶垫	Big instrument table silicone pad		pcs				✓		✓		✓		
91	2. 03. 202XZ003A1-1	大器械盘体	Big instrument table body		pcs				✓		✓		✓		
92	2. 03. 202XZ001A0-1	大器械盘前罩壳	Big instrument table front cover		pcs				✓		✓		✓		
93	2. 03. 202XZ004A0-1	手机挂架1	Handpiece holder 1	MC2, For 6-way syringe	pcs				✓		✓		✓		
94	2. 03. 202XZ005A0-1	手机挂架2	Handpiece holder 2	For fiber-optic handpiece/Satalec scaler	pcs				✓		✓		✓		
95	2. 03. 202XZ006A0-1	手机挂架3	Handpiece holder 3	Domestic turbine tubing /Low speed handpiece/For syringe	pcs				✓		✓		✓		
96	2. 03. 202XZ007A0-1	手机挂架4	Handpiece holder 4	Woodpecker,for EMS Scaler	pcs				✓		✓		✓		
97	1. 05. 150ZL080A2	强吸手柄总成(铝)	HV suction handle assembly (Aluminum)		pcs			✓	✓	✓	✓	✓	✓	✓	✓
98	1. 05. 150ZL090A1	弱吸手柄总成(铝)	Saliva ejector handle assembly (Aluminum)		pcs			✓	✓	✓	✓	✓	✓	✓	✓
99	3. 05. 350ZL084A1	强吸手柄开关小接头	HV suction handle switch small connector	Use for standard suction	pcs			✓	✓	✓	✓	✓	✓	✓	✓
100	3. 05. 350ZL094A1	弱吸手柄开关小接头	Saliva ejector handle switch small connector	Use for standard suction	pcs			✓	✓	✓	✓	✓	✓	✓	✓
101	3. 05. 350ZL093A1	弱吸手柄开关大接头	Saliva ejector handle switch big connector	Use for negative pressure suction	pcs			✓	✓	✓	✓	✓	✓	✓	✓
102	3. 05. 350ZL083A1	强吸手柄开关大接头	HV suction handle switch big connector	Use for negative pressure suction	pcs			✓	✓	✓	✓	✓	✓	✓	✓
103	2. 03. 201ZL418A0	弱吸管接头	Saliva ejector tubing connector		pcs			✓	✓	✓	✓	✓	✓	✓	✓
104	2. 03. 201ZL419A0	强吸管接头	HV suction tubing connector		pcs			✓	✓	✓	✓	✓	✓	✓	✓

S series

105	3.01.301CX075A0	强吸(开关)手柄接头	HV suction (switch) handle joint	Use for negative pressure suction	pcs		✓	✓	✓	✓	✓	✓
106	3.01.301CX076A0	弱吸(开关)手柄接头	Saliva ejector (switch) handle joint	Use for negative pressure suction	pcs		✓	✓	✓	✓	✓	✓
107	3.01.301CX007A0	强吸(开关)手柄	HV suction (switch) handle		pcs		✓	✓	✓	✓	✓	✓
108	3.01.301CX039A0	弱吸(开关)手柄	Saliva ejector (switch) handle		pcs		✓	✓	✓	✓	✓	✓
109	3.01.301CX042A0	弱吸接头(带O型圈)	Saliva ejector connector (with O-ring)		pcs		✓	✓	✓	✓	✓	✓
110	3.08.301AA011A0	弱吸转换接头	Saliva ejector adaptor		pcs		✓	✓	✓	✓	✓	✓
111	3.01.301CX008A1	强吸连通管	HV suction tubing	13.6*9.4*1700mm	pcs		✓	✓	✓	✓	✓	✓
112	3.01.301CX009A1	弱吸连通管	Saliva ejector tubing	10*6.0*1650mm	pcs		✓	✓	✓	✓	✓	✓
113	3.08.301AA033A0	灰色强吸管	Grey HV suction tubing	ID:17.5mm	meter		✓	✓	✓	✓	✓	✓
114	3.08.301AA034A0	灰色弱吸管	Grey saliva ejector tubing	ID:10.5mm	meter		✓	✓	✓	✓	✓	✓
115	2.03.201CX130A1	强弱吸过滤器盖	Suction filter cover		pcs		✓	✓	✓	✓	✓	✓
116	2.03.201CX136A2	强弱吸过滤器座	Suction filter holder		pcs		✓	✓	✓	✓	✓	✓
117	4.03.431AA15020	O型圈	O-ring (ID:15*2)	OD 15*2	pcs		✓	✓	✓	✓	✓	✓
118	4.03.431AA19180	O型圈	O-ring (OD:19*1.8)	ID 19*1.8	pcs		✓	✓	✓	✓	✓	✓
119	4.03.431AA20150	O型圈	O-ring (ID:20*1.5)	ID 20*1.5	pcs		✓	✓	✓	✓	✓	✓
120	4.03.431AA24150	O型圈	O-ring (ID:24*1.5)	ID 24*1.5	pcs		✓	✓	✓	✓	✓	✓
121	2.03.201CX137A0	强吸过滤网	HV Suction filter net		pcs		✓	✓	✓	✓	✓	✓
122	2.03.201CX138A0	弱吸过滤网	Saliva ejector filter net		pcs		✓	✓	✓	✓	✓	✓
123	3.08.301AA017A1	国产吸唾过滤器过滤网+杆	Domestic suction filter net + rod		pcs		✓	✓	✓	✓	✓	✓
124	3.08.301AA017A0	国产吸唾过滤器	Domestic suction filter		pcs		✓	✓	✓	✓	✓	✓
125	4.03.431AA15020	O型圈	O-ring (ID:15*2)	OD 15*2	pcs		✓	✓	✓	✓	✓	✓
126	4.03.431AA19180	O型圈	O-ring (OD:19*1.8)	ID 19*1.8	pcs		✓	✓	✓	✓	✓	✓
127	4.03.431AA35020	O型圈	O-ring (ID:35*2)	OD 35*2	pcs		✓	✓	✓	✓	✓	✓
128	2.03.201ZL417A0-1	治疗台小臂护盖	Instrument table connecting arm cover		pcs		✓	✓	✓	✓	✓	✓




















S series

129	2. 03. 201ZS401A3	助手位上盖	Assistant top cover		pcs		✓	✓	✓	✓	✓	✓
130	2. 03. 207CX220A3	口杯出水头	Cup filler		pcs		✓	✓	✓	✓	✓	✓
131	2. 03. 207CX221A2	痰盂出水头	Spittoon rinser		pcs		✓	✓	✓	✓	✓	✓
132	4. 03. 431AA10200	O型圈 (硅胶)	O-ring (ID:10*2)	OD10*2	pcs		✓	✓	✓	✓	✓	✓
133	2. 03. 201JT207A0	单向开关按钮	One-way switch button		pcs		✓	✓	✓	✓	✓	✓
134	2. 03. 201JT208A0	单向开关支架	One-way switch holder		pcs		✓	✓	✓	✓	✓	✓
135	2. 03. 201JT209A0	单向开关固定板	One-way switch fixed board		pcs		✓	✓	✓	✓	✓	✓
136	2. 03. 201JT211A0	脚踏叉簧	Foot control hookswitch		pcs		✓	✓	✓	✓	✓	✓
137	2. 03. 207JT001A7	U200脚开关上盖(蓝色丝印)	Foot control top cover(blue silkscreen)		pcs		✓	✓	✓	✓	✓	✓
138	2. 03. 254JT001A0	V2000脚开关上盖	Foot control top cover		pcs		✓	✓	✓	✓	✓	✓
139	2. 04. 251JT029A0	踏板脚开关提手	Lifting handle of pedal foot control		pcs		✓	✓	✓	✓	✓	✓
140	2. 06. 201JT206A1	复位弹簧	Reset spring		pcs		✓	✓	✓	✓	✓	✓
141	2. 03. 207JT002A8	脚开关底板	Foot control bottom plate		pcs		✓	✓	✓	✓	✓	✓
142	2. 03. 207JT009A2	小齿轮	Pinion		pcs		✓	✓	✓	✓	✓	✓
143	2. 03. 207JT021A1-1	拨杆组件	Lever assembly		pcs		✓	✓	✓	✓	✓	✓
144	2. 03. 207JT023A2	齿pcs	Rack		pcs		✓	✓	✓	✓	✓	✓
145	4. 02. 422QC05200	脚踏强磁	Foot pedal magnetic	φ5X2	pcs		✓	✓	✓	✓	✓	✓

S series

146	2. 06. 207JT025A4	脚开关扭簧一	Foot control torsional spring		pcs		✓	✓	✓	✓	✓	✓
147	2. 03. 207JT035A0-1	移动开关面板	Movement switch panel		pcs		✓	✓	✓	✓	✓	✓
148	2. 03. 251JT003A0	脚开关踏板	Foot control pedal		pcs		✓	✓	✓	✓	✓	✓
149	2. 06. 207JT036A0	脚开关拨杆限位块	Foot control lever limited block		pcs		✓	✓	✓	✓	✓	✓
150	4. 02. 422KG03006	微动开关	Micro switch		pcs		✓	✓	✓	✓		
151	4. 02. 422KG12003	微动开关	Micro switch		pcs		✓	✓	✓	✓	✓	✓
152	4. 02. 422KG14280	开关	Switch		pcs		✓	✓	✓	✓	✓	✓
153	1. 01. 121CX004F2	S注塑托盘支架总成	S Plastic injection tray support assembly	Flexible	pcs		✓		✓		✓	
154	2. 03. 251ZL161A2	V1000托盘	V1000 Tray		pcs		✓		✓		✓	
155	3. 05. 350ZL221A0	托盘硅胶垫	Tray silicone pad		pcs		✓		✓		✓	
156	1. 01. 111ZL011F0	下挂治疗台前臂连接杆	Instrument table holding arm (X-type)		pcs			✓		✓		✓
157	2. 04. 201SG116A5-1	上拉治疗台前臂连接杆	Instrument table holding arm (S-type)		pcs		✓		✓		✓	
158	4. 03. 430FC35500	含油轴承 (S)	Sleeve Bearing (S)		pcs		✓	✓	✓	✓	✓	✓
159	2. 04. 201ZL416A0-1	升降臂端盖	Lifting arm end cover		pcs		✓	✓	✓	✓	✓	✓
160	2. 04. 201ZS301A9	助手位组焊底板	Assistant bottom plate		pcs		✓	✓	✓	✓	✓	✓
161	2. 06. 201CX127A0	强吸下水大三通	HV Drainage 3-way link (big)		pcs		✓	✓	✓	✓	✓	✓
162	2. 06. 201XZ008A0	观片灯压片	Film viewer depressor		pcs		✓	✓	✓	✓	✓	✓

S series

163	2.06.203ZL120A3	电磁阀下挂固定板	X-type Solenoid valve fixer		pcs				✓		✓		✓
164	3.01.301CX025A0	子母扣头	Snap button (part 1)		pcs				✓	✓	✓	✓	✓
165	3.01.301CX026A0	子母扣座	Snap button (part 2)		pcs				✓	✓	✓	✓	✓
166	3.01.301CX034A3	小流量调节阀	Small flow control valve	φ3.2-φ2.3	pcs				✓	✓	✓	✓	✓
167	3.01.301CX065A0	小流量调节阀	Small flow control valve	2-φ3.2	pcs				✓	✓	✓	✓	✓
168	3.01.301CX044A2	分道阀	Separation valve	Two inlet&one outlet(φ3.2)	pcs				✓	✓	✓	✓	✓
169	3.01.301CX048A1	直角气锁开关	Instrument table air lock switch	φ3.2	pcs					✓		✓	
170	3.02.302CX004A1	气锁开关	Air lock switch	φ3.2	pcs				✓		✓		
171	4.03.430DQ18010	气锁开关螺母	Air lock switch nut	18*1	pcs				✓		✓		
172	3.01.301CX031A0	医生椅气弹簧	Air spring (Q7Q)		pcs				✓	✓	✓	✓	✓
173	2.03.201YY121A0-1	医生椅靠背罩壳	Backrest cover		pcs				✓	✓	✓	✓	✓
174	2.03.201YY105A1-1	医生座椅后罩A	Chair back cover A		pcs				✓	✓	✓	✓	✓
175	2.03.201YY106A1-1	医生座椅后罩B	Chair back cover B		pcs				✓	✓	✓	✓	✓
176	2.07.201YY112A6	靠背支架	Backrest holder		pcs				✓	✓	✓	✓	✓
177	2.07.201YY121A0	护士椅靠背	Assistant stool backrest		pcs				✓	✓	✓	✓	✓
178	3.01.301CX068A0	300铝合金五星脚	Aluminum alloy five-star feet		pcs				✓	✓	✓	✓	✓
179	3.01.301CX069A1	细头灰边PU脚轮	Small head grey PU caster		pcs				✓	✓	✓	✓	✓
180	3.01.301YY007A1	S系列护士椅底盘	S-series assistant chair base plate		pcs				✓	✓	✓	✓	✓
181	3.01.301YY008A0	S系列护士椅弯背	S-series assistant chair bent back		pcs				✓	✓	✓	✓	✓



















S series

182	3.01.301ZY068A0	西格医生椅底盘	SIGER dentist chair base plate		pcs		✓	✓	✓	✓	✓	✓
183	3.02.302CX005A1	四孔三联膜片阀	Triple-link membrane valve		pcs		✓	✓	✓	✓		
184	3.08.308CX005A1	四孔三联膜片阀膜片(老)	Triple-link membrane valve membrane (Old)		pcs		✓	✓	✓	✓		
185	3.04.307CX006A0	U痰盂装饰盖	U-series Spittoon decorative cover		pcs		✓	✓	✓	✓	✓	✓
186	4.02.422KG59520	船形开关	Rocket switch		pcs		✓	✓	✓	✓	✓	✓
187	4.04.432JF01700	压力表	Pressure gauge	GU-40	pcs		✓	✓	✓	✓	✓	✓
188	4.04.432JT00600	铜接头 (G1/2*5/16)	Copper connector (G1/2*5/16)	IBC8-04 (N03-C4025)	pcs		✓	✓	✓	✓	✓	✓
189	2.06.201ZL803A0	内窥镜支架	Intra-oral camera holder (inner bore 55mm)		pcs		✓	✓	✓	✓	✓	✓
190	3.01.301CX028A0	一级水过滤器 (不含接头)	Water filter grade 1 (without connector)		pcs		✓	✓	✓	✓	✓	✓
191	3.01.301CX028A1	一级水过滤器过滤网	Water filter net grade 1		pcs		✓	✓	✓	✓	✓	✓
192	3.01.301ZY006A0	万向头枕支架	Omni-directional headrest holder	Stainless Steel	pcs		✓	✓	✓	✓	✓	✓
193	3.01.301ZY005A0	头枕支架	Headrest holder	Stainless Steel	pcs		✓	✓	✓	✓	✓	✓
194	2.07.201ZY717A0	头枕堵盖	Headrest blanking cap		pcs		✓	✓	✓	✓	✓	✓
195	2.06.202XC008A5	大下挂器械盘PCB基板	Big X-type instrument table PCB base plate		pcs			✓		✓		✓
196	3.04.307JT015A1	S气控脚开关阀	Pneumatic foot control valve	(lever type)	pcs		✓	✓	✓	✓		
197	3.05.354JT401A0	V2000踏板脚开关调速阀	Foot pedal switch speed control valve (press type)	(pedal type)	pcs		✓	✓	✓	✓		
198	3.02.302ZC008A1	S60侧箱水气分配器总成	Water & air distribution assembly in hydro group		pcs		✓	✓	✓	✓	✓	✓

S series

199	3. 08. 301AA206A0	METASYS选位阀	METASYS selection valve		pcs		✓	✓	✓	✓	✓	✓
200	3. 08. 301AA003A0	污物分离器	Contaminant separator		pcs		✓	✓	✓	✓	✓	✓
201	3. 08. 301AA004A0	干吸截止阀	Dry suction stop valve		pcs		✓	✓	✓	✓	✓	✓
202	3. 08. 301AG100A0	黑色光导棒(啄木鸟)	Black optical rod (Woodpecker)		pcs		✓	✓	✓	✓	✓	✓
203	3. 08. 308AE009A0	啄木鸟光固化遮光片	Woodpecker curing light brake sack		pcs		✓	✓	✓	✓	✓	✓
204	4. 05. 401SG025F4	脚踏线	Foot control cable	10 core L=2000	pcs		✓	✓	✓	✓		
205	4. 05. 401SG026F5	S90脚踏线	Foot control cable		pcs						✓	✓
206	1. 01. 111ZC004A0	热水器总成	Water heater assembly		pcs		✓	✓	✓	✓	✓	✓
207	1. 01. 121ZY003F0	90头枕固定座总成	Headrest fixer assembly		pcs		✓	✓	✓	✓	✓	✓
208	2. 07. 203ZY805A4	左扶手	Left armrest		pcs		✓	✓	✓	✓	✓	✓
209	1. 08. 191ZY001F1	S活动扶手	S Mobile armrest	With 30 armrest press block	pcs		✓	✓	✓	✓	✓	✓
210	1. 08. 191ZY003A0	S活动左扶手	S Mobile left armrest	With 30 armrest press block	pcs		✓	✓	✓	✓	✓	✓
211	1. 05. 150ZY005F0	U500上平行盖总成	Chair upper cover assembly		pcs		✓	✓	✓	✓	✓	✓
212	2. 03. 250ZY808A3	U500底座盖板	Base cover board		pcs		✓	✓	✓	✓	✓	✓
213	2. 03. 250ZY803A3	U500下平行盖	Chair lower cover		pcs		✓	✓	✓	✓	✓	✓
214	2. 03. 201CX052A1	PCB小盒(底壳+上盖)	Small PCB box(bottom base + top cover)		pcs		✓	✓	✓	✓	✓	✓
215	2. 03. 201CX053A1	PCB大盒(底壳+上盖)	Big PCB box(bottom base + top cover)		pcs		✓	✓	✓	✓	✓	✓
216	2. 03. 201ZY710A0-1	头枕锁紧手柄	Headrest locker		pcs		✓	✓	✓	✓	✓	✓

S series

217	2. 06. 250CX601A0	过滤器转换接头	Filter adaptor		pcs			✓	✓	✓	✓	✓	✓
218	2. 04. 250ZY430A3	U500铝靠背	U500 Aluminium backrest		pcs			✓	✓	✓	✓	✓	✓
219	2. 07. 250ZY432A1	U500靠背卡扣盖	U500 backrest buckle cover		pcs			✓	✓	✓	✓	✓	✓
220	2. 07. 250ZY435A0	加长靠背子母扣	Extended backrest Snap button		pcs			✓	✓	✓	✓	✓	✓
221	2. 20. 220DL008A1	口腔灯把手固定簧卡	Operating light handle spring fixer	φ1.2mm	pcs			✓	✓	✓	✓	✓	✓
222	3. 01. 301CX015A1	带回气气开关	Air switch with air release function	φ3.2	pcs			✓	✓	✓	✓	✓	✓
223	3. 01. 301CX018A1	弯头三用枪	Angled syringe	(φ3.2)	pcs			✓	✓	✓	✓	✓	✓
224	3. 01. 301CX027A0	强吸嘴	HV suction nozzle		pcs			✓	✓	✓	✓	✓	✓
225	3. 01. 301CX029A1	直头三用枪	Straight syringe	φ3.2	pcs			✓	✓	✓	✓	✓	✓
226	3. 01. 301CX051A1	管道节流阀	Throttling valve	Copper	pcs			✓	✓	✓	✓	✓	✓
227	3. 01. 301CX062A0	下水组合	Launching assembly		pcs			✓	✓	✓	✓	✓	✓
228	3. 01. 301WG003A1	弱吸发生器	Saliva ejector generator	Copper	pcs			✓	✓	✓	✓	✓	✓
229	3. 01. 301WG004A1	强吸发生器	HV suction generator	Copper	pcs			✓	✓	✓	✓	✓	✓
230	3. 02. 302CX002A0	总水气机械阀	Water & air mechanical valve	JCF-1V	pcs			✓	✓	✓	✓	✓	✓
231	3. 02. 302ZC011A1	60三通总成	60 3-way link assembly	3*φ3.2	pcs			✓	✓	✓	✓	✓	✓
232	3. 02. 302ZC012A0	4.8三通总成	3-way link assembly φ4.8	3*φ4.8	pcs			✓	✓	✓	✓	✓	✓
233	3. 02. 303ZC003A1	30五通-三通组合总成	5-links & Triple-links assembly	1*φ4.8+6*φ3.2	pcs			✓	✓	✓	✓	✓	✓
234	3. 04. 303XP001A1	气控电阀	Pneumatic electric valve	φ3.2	pcs			✓	✓	✓	✓	✓	✓


S series

235	3.02.302CX018A0	膜片式可调压力开关	Membrane adjustable pressure switch		pcs		✓	✓	✓	✓		
236	3.04.307CX004A2	气弹簧 (Q7Q)	Air spring (Q7Q)	485*1500N	pcs		✓	✓	✓	✓	✓	✓
237	3.04.307ZL012A4	U200治疗台分配器	U200 Instrument table distributor		pcs		✓	✓	✓	✓	✓	✓
238	3.04.309XP003A1	六角压水阀(进口膜片)	Hexangular water control valve (imported diaphragm)	φ3.2	pcs		✓	✓	✓	✓		
239	3.04.310CX004A0	水气快接头	Water & air quick connector		pcs		✓	✓	✓	✓	✓	✓
240	3.08.301AA120A0	国产选位阀(配控制电位器)	Domestic selection valve(With control potentiometer)		pcs		✓	✓	✓	✓	✓	✓
241	3.08.301XP013A0	三用枪喷嘴	Syringe nozzle		pcs		✓	✓	✓	✓	✓	✓
242	4.01.416BA10521	保险丝	Fuse	T10A/250V 5x20	pcs		✓	✓	✓	✓	✓	✓
243	4.01.416JX11624	继电器(负压用)	Relay(negative pressure)		pcs		✓	✓	✓	✓	✓	✓
244	4.02.422KG10003	微动开关 (三pcs脚)	Micro switch (3 pins)		pcs		✓	✓	✓	✓	✓	✓
245	4.02.422KG10004	微动开关(无弹片)	Micro switch (Without shrapnel)		pcs		✓	✓	✓	✓	✓	✓
246	2.03.250ZY801A2	U500座椅前罩	Chair front cover		pcs		✓	✓	✓	✓	✓	✓
247	2.02.250CX305A1	U500侧箱座椅连接盖	Connecting cover between hydro group & chair		pcs		✓	✓	✓	✓	✓	✓
248	2.02.250ZY804A2	U500座椅后盖	Chair back cover		pcs		✓	✓	✓	✓	✓	✓
249	2.02.250ZY805A2	座垫前护壳	Cushion front cover		pcs		✓	✓	✓	✓	✓	✓
250	2.03.250ZY702A2	导轨滑块	Track sliding block		pcs		✓	✓	✓	✓	✓	✓
251	6.04.650ZY701A1	导轨	Track		pcs		✓	✓	✓	✓	✓	✓
252	2.07.250ZY414A6	摆杆	Oscillating bar		pcs		✓	✓	✓	✓	✓	✓

S series

253	2. 04. 250ZY405A0	U500弯板	Curved plate		pcs		✓	✓	✓	✓	✓	✓
254	4. 04. 432ZL60600	直通快接 φ6	Straight connector φ6	φ6-φ6	pcs		✓	✓	✓	✓	✓	✓
255	1. 01. 111CX002F1	侧箱门板总成(瓷白)	Hydro group cover assembly (white)		pcs		✓	✓	✓	✓	✓	✓
256	2. 03. 201CX404A0-1	侧箱活动盖	Hydro group mobile cover housing		pcs		✓	✓	✓	✓	✓	✓
257	2. 03. 201CX403A0-1	侧箱固定盖	Hydro group fixed cover housing		pcs		✓	✓	✓	✓	✓	✓
258	1. 01. 111CX403F1	侧箱固定盖(瓷白加工)	Hydro group fixed cover housing (White)		pcs		✓	✓	✓	✓	✓	✓
259	1. 01. 111ZC003A0	气缸总成	Air cylinder assembly		pcs				✓	✓	✓	✓
260	1. 01. 111ZC006A1	漱口盆总成	Spittoon assembly		pcs		✓	✓	✓	✓	✓	✓
261	1. 01. 111ZC020F3	助手位调节臂总成(瓷白)	Assistant adjusting arm assembly(white)		pcs				✓	✓	✓	✓
262	2. 03. 201ZS402A0-1	调节臂后盖	Adjusting arm back cover		pcs				✓	✓	✓	✓
263	2. 03. 201ZS403A0-1	调节臂前盖	Adjusting arm front cover		pcs				✓	✓	✓	✓
264	2. 03. 201CX205A6	助手位转盘	Assistant rotating housing		pcs		✓	✓	✓	✓	✓	✓
265	2. 04. 201ZL101B0	治疗台转臂	Instrument table rotating arm		pcs		✓	✓	✓	✓	✓	✓
266	3. 01. 301CX004A1	水源转换开关	Water source switch	φ3.2	pcs		✓	✓	✓	✓	✓	✓
267	3. 01. 301CX061A2	8寸钢珠滑轨	8" steel ball sliding rail		付		✓	✓	✓	✓	✓	✓
268	3. 01. 301CX311A0	S蒸馏水瓶	Distilled water bottle		pcs		✓	✓	✓	✓	✓	✓
269	1. 01. 111ZC029A0	RT蒸馏水瓶盖总成	RT Distilled water bottle cap assembly		pcs		✓	✓	✓	✓	✓	✓




















S series

270	3. 01. 301CX019A0	痰盂盆	Spittoon		pcs		✓	✓	✓	✓	✓	✓
271	2. 03. 201CX405A0	痰盂漏斗	Spittoon funnel		pcs		✓	✓	✓	✓	✓	✓
272	2. 03. 201CX407A3-1	底箱罩壳	Base box cover housing		pcs		✓	✓	✓	✓	✓	✓
273	2. 04. 220DD009A0-1	S侧箱大臂孔堵盖	Hydro group arm hole blanking cap		pcs		✓	✓	✓	✓	✓	✓
274	4. 02. 423CB230F3	S30/60/90 230V变压器	Transformer(230V)	230V	pcs		✓	✓	✓	✓	✓	✓
275	4. 02. 423CB120A0	S30/60/90 120V变压器	Transformer(120V)	120V	pcs		✓	✓	✓	✓	✓	✓
276	2. 04. 201ZL117A8	手术灯杆	Light pole	L=915	pcs		✓	✓	✓	✓	✓	✓
277	2. 04. 201ZL120B2	内窥镜灯杆	Oral camera lamp pole	L=915	pcs		✓	✓	✓	✓	✓	✓
278	1. 02. 117ZC008A0	U200侧箱电磁阀总成 (含接头和支架)	Solenoid valve assembly in hydro group (With connector and holder)		pcs		✓	✓				
279	1. 07. 117ZC008A0	U200侧箱电磁阀总成 (含接头)	Solenoid valve assembly in hydro group (With connector)		pcs		✓	✓				
280	3. 04. 307CX001A0	中型电磁阀 (不含接头)	Medium-size solenoid valve (Without connector)		pcs		✓	✓				
281	1. 05. 150CX007A0	U500微型电磁阀总成	Micro solenoid valve assembly	φ3.2	pcs		✓	✓				
282	1. 05. 150CX008A0	U500三联电控阀总成	Triple-link electrical control valve assembly	φ3.2	pcs		✓	✓	✓	✓		
283	1. 07. 112ZC010A0	膜片阀用电磁阀总成 (维修用)	Solenoid valve assembly for membrane valve (for service)		pcs		✓	✓	✓	✓		
284	1. 01. 112ZC015A0	60器械盘小流量电磁阀总成	Instrument table small flow solenoid valve assembly		pcs				✓	✓		
285	1. 01. 111ZC025A0	侧箱小流量电磁阀总成	Hydro group small flow solenoid valve assembly		pcs				✓	✓	✓	✓
286	1. 01. 111ZC016A1	侧箱电磁阀总成 (冲盂用)	Solenoid valve assembly for hydro group (For rinsing)		pcs				✓	✓	✓	✓

S series

287	1. 01. 111ZC013A1	气锁电磁阀总成	Airlock solenoid valve assembly	φ3.2	pcs									✓	✓
288	1. 01. 111ZC014A1	比例电磁阀总成	Proportional solenoid valve assembly		pcs									✓	✓
289	1. 01. 111ZC026A0	小电磁阀总成	Small solenoid valve assembly	φ3.2	pcs									✓	✓
290	3. 01. 301CX001A3	S90器械盘组合阀	S90 group valve assembly for instrument table		pcs									✓	✓
291	4. 05. 401SG122F1	负压泵电源线	Suction pump power line	2*0.75 L=2300	pcs		✓	✓	✓	✓	✓	✓	✓	✓	✓
292	1. 01. 111CX2094A	出水口座(加工)	Water outlet base		pcs		✓	✓	✓	✓	✓	✓	✓	✓	✓
293	4. 03. 431AA03612	O型圈	O-ring	ID 3.6*1.2	pcs		✓	✓	✓	✓	✓	✓	✓	✓	✓
294	4. 03. 431AC31542	出水口座垫圈(硅胶)	Water outlet base gasket (silica gel)	31.5*40*2	pcs		✓	✓	✓	✓	✓	✓	✓	✓	✓
295	2. 04. 201CX209A1	出水头	Water outlet		pcs		✓	✓	✓	✓	✓	✓	✓	✓	✓
296	2. 03. 201ZY714A3	头枕罩壳	Armrest cover		pcs		✓	✓	✓	✓	✓	✓	✓	✓	✓
297	4. 03. 430DQ08700	膨胀螺丝	Setscrew	8*70	pcs		✓	✓	✓	✓	✓	✓	✓	✓	✓
298	3. 04. 307CX003A0	φ50软管螺母(上连接套)	φ50 Soft tube nut(up adapter sleeve)		pcs		✓	✓	✓	✓	✓	✓	✓	✓	✓
299	2. 02. 230YD015A0	灯杆装饰盖	Operating light pole decorative cover		pcs		✓	✓	✓	✓	✓	✓	✓	✓	✓
300	4. 03. 432AA13190-01	喉箍	Hose clamp	13-19mm	pcs		✓	✓	✓	✓	✓	✓	✓	✓	✓
301	4. 03. 432AA28310	塑料卡箍	Plastic clamp	28.8-31.1/For ID 25 tube use	pcs		✓	✓	✓	✓	✓	✓	✓	✓	✓
302	4. 03. 432AA18200	塑料卡箍	Plastic clamp	18.8-20.6/For ID 14 tube use	pcs		✓	✓	✓	✓	✓	✓	✓	✓	✓
303	4. 03. 432AA20210	塑料卡箍	Plastic clamp	20.1-21.9/For ID 17 tube use	pcs		✓	✓	✓	✓	✓	✓	✓	✓	✓












S series

304	4. 03. 432AA17190	塑料卡箍	Plastic clamp	17.2-19.2/For HV suction、saliva ejector generator use	pcs		✓	✓	✓	✓	✓	✓
305	2. 01. 250ZY002A0	手术椅按键连接板	Chair button connection board	SSYAJ-CON-A0	pcs		✓	✓	✓	✓	✓	✓
306	2. 01. 250ZY001A1	手术椅按键板	Chair button control board	SSYAJ-01-A1	pcs		✓	✓	✓	✓	✓	✓
307	3. 05. 350MB007A0	座椅控制面板	Chair control board		pcs		✓	✓	✓	✓	✓	✓
308	2. 03. 252ZY411A0	座椅控制器上盖	Chair control board top cover		pcs		✓	✓	✓	✓	✓	✓
309	4. 04. 432JF00803	气调压过滤器(配Y40压力表无支架)	Air filter with pressure adjustment (with Y40 pressure meter, without holder)	AFR2000M1,ZH-027C	pcs		✓	✓	✓	✓	✓	✓
310	2. 03. 220DL011A0	LED口腔灯透明罩	LED operating light transparent cover		pcs		✓	✓	✓	✓	✓	✓
311	2. 03. 220DL012A0	LED口腔灯前罩	LED operating light front cover		pcs		✓	✓	✓	✓	✓	✓
312	2. 03. 220DL018A0	LED口腔灯后罩	LED operating light back cover		pcs		✓	✓	✓	✓	✓	✓
313	2. 03. 220DL016A0	LED口腔灯后罩卡扣壳	LED operating light back cover buckle shell		pcs		✓	✓	✓	✓	✓	✓
314	2. 03. 220DL017A0	LED口腔灯后罩卡扣盖	LED operating light back cover buckle cover		pcs		✓	✓	✓	✓	✓	✓
315	2. 20. 220DL025A1	口腔灯后盖固定簧卡	LED operating light back cover spring fixer	Φ1.0mm	pcs		✓	✓	✓	✓	✓	✓
316	2. 03. 220DL019A1	LED口腔灯接头悬臂	LED operating light connector cantilever		pcs		✓	✓	✓	✓	✓	✓
317	4. 03. 431AA25021	O型圈O-RING NBR 70 BLACK NORMAL	O-RING NBR 70 BLACK NORMAL	Inner∅ 25*2	pcs		✓	✓	✓	✓	✓	✓
318	2. 03. 220DL020A1	LED口腔灯上接头	LED operating light up connector		pcs		✓	✓	✓	✓	✓	✓
319	2. 20. 220DL024A2	上接头转轴(佛山)	Up connector rotary shaft (Foshan)		pcs		✓	✓	✓	✓	✓	✓
320	4. 03. 431AA27250	O型圈	O-ring	Inner∅ 27*2.5	pcs		✓	✓	✓	✓	✓	✓
321	2. 03. 220DL026A1	口腔灯装饰圈	Operating light decorative ring		pcs		✓	✓	✓	✓	✓	✓
322	2. 03. 222DS003A0	反光镜片	Reflector		pcs		✓	✓	✓	✓	✓	✓

S series

323	1. 07. 120DL022A1	LED口腔灯把手总成	LED operating light handle assembly		pcs		✓	✓	✓	✓	✓	✓
324	1. 07. 120DL014A2	LED悬臂转轴总成	LED cantilever rotary shaft assembly		pcs		✓	✓	✓	✓	✓	✓
325	2. 03. 220DL013A0	LED口腔灯光阑	LED operating light lens frame		pcs		✓	✓	✓	✓	✓	✓
326	2. 03. 220DL001A0	LED口腔灯聚光镜片	LED operating light focusing lens		pcs		✓	✓	✓	✓	✓	✓
327	2. 20. 220DL005A1	LED口腔灯散热片	LED operating light cooling fin		pcs		✓	✓	✓	✓	✓	✓
328	2. 20. 220DL002A2	LED口腔灯镜片遮光膜	LED operating light lens anti-glare membrane		pcs		✓	✓	✓	✓	✓	✓
329	2. 20. 221DP004A0	把手转轴	Handle rotary shaft		pcs		✓	✓	✓	✓	✓	✓
330	4. 04. 432JF01600	三通	Triple-way link		pcs		✓	✓	✓	✓	✓	✓
331	4. 04. 432JF03960	隔膜止回阀	Membrane check valve		pcs		✓	✓	✓	✓	✓	✓
332	4. 04. 432JF012A1	简易减压阀	Simple pressure reducing valve	3,2	pcs		✓	✓	✓	✓	✓	✓
333	4. 04. 432JF01800	内螺纹直通	Internal thread straight connection		pcs		✓	✓	✓	✓	✓	✓
334	4. 05. 401SG044F2	电马达电源线	Electric motor power line	L=420	pcs		✓	✓	✓	✓	✓	✓
335	7. 701ZY410A2	PVC透明脚垫	PVC transparent foot pad	740*350	pcs		✓	✓	✓	✓	✓	✓
336	4. 05. 440TT04400	铜管套 φ4.4	Copper tube cap φ4.4	φ4.4	pcs		✓	✓	✓	✓	✓	✓
337	4. 05. 440TT04800	铜管套 φ4.8	Copper tube cap φ4.8	φ4.8	pcs		✓	✓	✓	✓	✓	✓
338	4. 05. 440TT06400	铜管套 φ6.4	Copper tube cap φ6.4	φ6.4	pcs		✓	✓	✓	✓	✓	✓
339	4. 04. 432JD32321	对接头	Butt joint	φ3.2-φ3.2	pcs		✓	✓	✓	✓	✓	✓

S series

340	4. 04. 432JD48481	对接头	Butt joint	φ4.8-φ4.8	pcs			✓	✓	✓	✓	✓	✓
341	4. 04. 432JZ26500	铜接嘴	Copper joint	φ2.6-M5	pcs			✓	✓	✓	✓	✓	✓
342	4. 04. 432JZ32500	铜接嘴	Copper joint	φ3.2-M5	pcs			✓	✓	✓	✓	✓	✓
343	4. 04. 432JZ37500	铜接嘴	Copper joint	φ3.7-M5	pcs			✓	✓	✓	✓	✓	✓
344	4. 04. 432JZ48500	铜接嘴	Copper joint	φ4.8-M5	pcs			✓	✓	✓	✓	✓	✓
345	4. 05. 440PV01400	耐寒管	Cold-resistant tube	Inner φ 14 5 meters per roll	meter			✓	✓	✓	✓	✓	✓
346	3. 04. 307CX002A1	φ 50柔性软管	Soft tube φ50	φ50×60MM×1000MM	pcs			✓	✓	✓	✓	✓	✓
347	4. 05. 440PV02500	PVC加筋管	PVC rifflid tube	Inner 25 (L=1.78 meter)	pcs			✓	✓	✓	✓	✓	✓
348	4. 05. 440TG01600	柔性软管	Soft tube	φ16	meter			✓	✓	✓	✓	✓	✓
349	4. 05. 440PA47551	亚大PA管	PA Tube	4*0.75*500(White)	meter			✓	✓	✓	✓	✓	✓
350	4. 05. 440PA47553	亚大PA管	PA Tube	4*0.75*500(Yellow)	meter			✓	✓	✓	✓	✓	✓
351	4. 05. 440PA47552	亚大PA管	PA Tube	4*0.75*500(Blue)	meter			✓	✓	✓	✓	✓	✓
352	4. 05. 440PA64031	亚大PA管	PA Tube	6*4*300(White)	meter			✓	✓	✓	✓	✓	✓
353	4. 05. 440PA64032	亚大PA管	PA Tube	6*4*300(Blue)	meter			✓	✓	✓	✓	✓	✓
354	4. 05. 440PA86021	亚大PA管	PA Tube	8*6*200(White)	meter			✓	✓	✓	✓	✓	✓
355	4. 05. 440PA86022	亚大PA管	PA Tube	8*6*200(Blue)	meter			✓	✓	✓	✓	✓	✓